

INTRODUCING THE NEXT GENERATION OF ALL CERAMICS



Dentistry by
Dr. Michele Chalmers

**PREPAY 10 UNITS
AT \$109. SAVE \$300
(\$30 PER UNIT)**

Terms and conditions: Only applies to one unit restorations and not bridges. Shipping and tax not included. Cannot be combined with any other offer or discounts. Program expires Dec. 31, 2019.

Modern Dental's IPS e.max® ZirCAD® Prime delivers premium esthetics and exceptional strength

As the next generation of monolithic ceramics, IPS e.max ZirCAD Prime delivers an unprecedented combination of esthetics and strength. Zirconia durability and innovative gradient technology which provides seamless shade progression, make it an ideal material for single crowns in the smile zone and long-span bridges.

IPS e.max® ZirCAD® Prime	
Cementation: Conventionally cemented	Material: Oxide ceramic
Indications: Single crowns, to 14-unit bridges	Flexural strength: 1200 MPa
Preparations: 0.8 mm for anteriors and 1.0 mm for posteriors	Esthetics: Comparable to glass-ceramics

Experience a new level of all ceramics.

Contact Modern Dental at 810.534.9273 for more information and to schedule a case.