



# Guide to Cementation Selection





## Variolink® II


ESTHETIC RESIN CEMENT

### Features

- Proven long term clinical success
- Superior strength with outstanding esthetics
- (6) Light cure base shades and (2) catalyst shades
- (3) Viscosities
- Highest Radiopacity in its category



Variolink® II's excellent radiopacity permits the clear distinction of the restoration, cement, and caries on x-rays.




## Appeal™

ESTHETIC RESIN CEMENT

### Features

- Unique chemistry provides fast and easy clean up
- Highly translucent "Value" shades enhance the esthetics of the final restoration
- (7) Light Cure, (3) Dual cure values
- Microfill chemistry provides low wear rates and high polishability
- Amine free for outstanding color stability





## Multilink®

UNIVERSAL ADHESIVE CEMENT

### Features

- High bond strengths provide the ultimate security for all types of restorations.
- Universal applications make Multilink an ideal cement for indirect restorations (inlays, onlays, crowns, bridges, and posts) made of:
  - Metal and metal ceramics (PFM)
  - Metal-free (all-ceramics)
  - Strengthened core ceramics (Aluminum oxide, zirconia)
  - Composites
- Multilink applies quickly and cleans up easily.
- Self-etching and self-curing
- Three shades with distinct translucencies to meet all clinical indications.



## Vivaglass® CEM

GLASS IONOMER CEMENT

### The Easiest Method of Conventional Cementation

#### Features


- Convenient injection capsule
- Predosed mixing ratio
- Flexible applicator tip
- Low film thickness
- High radiopacity
- Continuous fluoride release
- Ideal capsule volume - enough for 1-2 crowns



# Selecting the Right Resin Cement

Recent developments in restorative materials and cements continue to support practitioners' ability to deliver aesthetic dental care. The proliferation of these therapeutic alternatives can be challenging, particularly when the clinician is

treatment planning and selecting the proper techniques, ceramic systems, and cements. The following guide highlights the necessary steps and benefits associated with three cements that are currently available at Ivoclar Vivadent.



When the clinician has three quality cements from which to choose, as shown herein, the ideal selection for a given indication often amounts simply to a matter of personal preference. Multilink, Appeal, and Variolink II can all be used chairside for optimal aesthetic results.

(Photographs and dentistry courtesy of Chris Ramsey, DDS.)  
 \*Cast gold crown was also seated with Multilink.

## Appeal™ ESTHETIC RESIN CEMENT

### Appeal Light Cure Anterior Assortment – 574341

#### Appeal Light Cure Base Values - 1x1.5 ml

574340 High Value +3 (Bleach)  
574348 High Value +2 (B1)  
574347 High Value +1 (Transparent)  
574346 Medium Value 0 (Clear)  
574344 Low Value -1 (A2)  
574345 Low Value -2 (A3)  
574342 Low Value -3 (A4)

#### Appeal Base Try-In Values - 1x1.5 ml

574332 High Value +3  
574339 High Value +2  
574338 High Value +1  
574337 Medium Value 0  
574334 Low Value -1  
574336 Low Value -2  
574333 Low Value -3

1 Monobond® S  
1 Liquid Strip®  
1 Total Etch®  
1 Excite® Bottle  
1 Excite® Applicators  
1 Brush holder  
Various canulas and brushes  
1 Vivapad  
1 Mixing pad  
1 Appeal Light Cure Shade Guide  
Instructions for use and flowchart



### Appeal Dual Cure Posterior Assortment - 574343

#### Appeal Dual Cure Values – 1x3 ml

574330 High Value +1 (Transparent)  
574331 Medium Value 0 (A1)  
574329 Low Value -1 (A4)

1 Monobond® S  
1 Liquid Strip®  
1 Total Etch®  
1 Excite® Bottle  
1 Excite® Applicators  
1 Brush holder  
Various canulas and brushes  
1 Vivapad  
1 Mixing pad  
1 Appeal Dual Cure Shade Guide  
Instructions for use and flowchart



## Variolink® II ESTHETIC RESIN CEMENT

### Esthetic Cementation System - 558946

#### Base – 1x2.5 gm

558949 Bleach XL  
558950 Transparent  
558951 White(110/A1)  
558954 White Opaque  
558952 Yellow(210/A3)  
558953 Brown(340/A4)



#### Catalyst – 1x2.5 gm

558955 High Viscosity Yellow(210/A3)  
558956 Low Viscosity Yellow(210/A3)  
558957 High Viscosity Transparent  
558958 Low Viscosity Transparent  
559224 Precision Applicator Tips

1 Total Etch®, 2 gm  
1 Bottle Excite®, 5 gm  
1 Monobond S®, 5 gm  
1 Liquid Strip, 2.5 gm  
5 Variolink®II Applicator Tips  
1 Brush Holder  
1 Vivapad®  
1 Mixing Pad  
2 VivaSticks®  
50 Disposable Brushes  
Various Cannulas

### Variolink II Try-In Standard Package - 558959

6 Syringes Variolink II Try-In; Shades Bleach XL(010), Transparent, White(110/A1), Yellow(210/A3), Brown(340/A4), White Opaque

#### Try-In – 1x2.5 gm

558960 Bleach XL(010)  
548793 Transparent  
548794 White(110/A1)  
548795 White Opaque  
548796 Yellow(210/A3)  
548797 Brown(340/A4)

### Intro Package - 560464

1 Syringe Variolink II Base 2.5 gm, Yellow(210/A3)  
1 Syringe Variolink II Catalyst 2.5 gm; low viscosity; Transparent  
1 Total Etch  
5 Excite, 0.1 gm  
1 Monobond S, 5 gm  
1 Liquid Strip, 2.5 gm  
Various Accessories

### Assortment - 558947

6 Syringes Variolink II Base 2.5 gm in shades Bleach XL(010), Transparent, White(110/A1), Yellow(210/A3), Brown(340/A4), White Opaque  
2 Syringes Variolink II Catalyst 2.5 gm; high/low viscosity; Transparent  
2 Syringes Variolink II Catalyst 2.5 gm; high/low viscosity; Yellow(210/A3)  
Various Accessories

### Variolink Ultra Standard Package - 538319

1 Syringe Ultra Base, 3 gm, shade Yellow(210/A3)  
1 Syringe Ultra Catalyst, 3 gm, shade Yellow(210/A3)

## Multilink® UNIVERSAL ADHESIVE CEMENT

### Multilink System Pack

576862 Yellow Shade  
576869 Transparent Shade  
576864 White Opaque Shade



#### Multilink System Pack Contains:

(1) Multilink Double Syringe, 6 gr  
(1) Multilink Primer A, 3 gr  
(1) Multilink Primer B, 3 gr  
(1) Multilink Metal primer, 5 gr  
Accessories - (50)brushes, (1)mixing pad, instructions for use

576825 Multilink Primer A, Primer B, contains 3 gr bottle  
576827 Multilink Metal Primer, 5 gr

#### Multilink Refills - 1x6 gr

576865 Yellow Shade  
576860 Transparent Shade  
576866 White Opaque Shade



## Vivaglass® CEM GLASS IONOMER CEMENT

544414 Standard Pack  
25 Capsules Vivaglass CEM

544413 Capsule Injector



USA  
175 Pineview Drive, Amherst, NY 14228  
716.691.0010 • 1.800.533.6825  
716.691.2285 Fax  
[www.ivoclarvivadent.us.com](http://www.ivoclarvivadent.us.com)

[www.ivoclarvivadent.us.com](http://www.ivoclarvivadent.us.com)

Call us toll free at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.  
© 2004 Ivoclar Vivadent, Inc. Appeal, Multilink, Variolink and Vivaglass are trademarks of Ivoclar Vivadent, Inc.

CANADA  
2785 Skymark Ave. Unit 1, Mississauga, Ontario L4W 4Y3  
905.238.5700 • 1.800.263.8182 (Inside Canada)  
905.238.5711 Fax  
[www.ivoclarvivadent.ca](http://www.ivoclarvivadent.ca)

100% CUSTOMER SATISFACTION  
GUARANTEED!

ivoclar  
vivadent  
clinical