


Provide Flexible and Durable Solutions for Your Patients!

ORTHO COLOR CHART



ICONS



MOUTHGUARD COLOR CHART



SPORTGUARD SPECIFICATION

Type	Thickness	Application	Remarks
Junior	3mm	Non-contact sports (junior)	Gymnastics, Figure Skating, Skateboard
Standard	5mm	Non-contact sports (adult)	Intense Cycling, Rock Climbing, Weight-lifting
Senior	6mm	Sports with some defree of contact	Wakeboarding, Baseball, Squash
Elite	6mm + 3mm at anterior adjusted to correct occlusion	Contact sports where objects with a large surface area may be received by the jaw	Soccer, Basketball, Wrestling
Professional	6mm + 3mm hard labial insert	Heavy contact sport where hard, pointed blows may be recieved by the jaw	Rugby, Boxing, Hockey